

BCL2L13 Antibody

Catalog No: #32164

Package Size: #32164-1 50ul #32164-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BCL2L13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total BCL2L13 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human BCL2L13 (NP_056182.2).
Target Name	BCL2L13
Other Names	BCL2L13;BCL-RAMBO;Bcl2-L-13;MIL1
Accession No.	Uniprot:Q9B XK5GeneID:23786
Uniprot	Q9B XK5
GeneID	23786
SDS-PAGE MW	48kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using BCL2L13 antibody.

Background

This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene.

Note: This product is for in vitro research use only