

Osteocrin Antibody HRP Conjugated

Catalog No: #C08057H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Osteocrin Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | IHC-P IHC-F ICC |
| Species Reactivity | Hu Ms Rt |
| Immunogen Description | KLH conjugated synthetic peptide aa 83-133 133 derived from human Osteocrin |
| Conjugates | HRP |
| Target Name | Osteocrin |
| Other Names | MUSCLIN; Osteocrin; OSTN |
| Accession No. | Swiss-Prot#P61366NCBI Gene ID344901 |
| Uniprot | P61366 |
| GeneID | 344901; |
| Excitation Emission | N A |
| Cell Localization | Secreted |
| Concentration | 1mg ml |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Appears to modulate osteoblastic differentiation. Could also function as an autocrine and paracrine factor linked to glucose metabolism in skeletal muscle (By similarity).

Note: This product is for in vitro research use only