

HSP27(Ab-82) Antibody

Catalog No: #21240

Package Size: #21240-1 50ul #21240-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HSP27(Ab-82) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total HSP27 protein.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around aa. 80~84 (Q-L-S-S-G) derived from Human HSP27.
Target Name	HSP27
Other Names	28 kDa heat shock protein; Estrogen-regulated 24 kDa protein; Growth-related 25 kDa protein; HSP 27; HSP25
Accession No.	Swiss-Prot: P04792NCBI Protein: NP_001531.1
Uniprot	P04792
GeneID	3315;
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 27kd


Western blotting: 1:500~1:1000

Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from K562 and HeLa cells using HSP27(Ab-82) Antibody #21240.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using HSP27(Ab-82) Antibody #21240(left) or the same antibody preincubated with blocking peptide(right).

Background

Involved in stress resistance and actin organization.

Note: This product is for in vitro research use only