

APOH Antibody

Catalog No: #32237

Package Size: #32237-1 50ul #32237-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	APOH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total APOH protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human APOH (NP_000033.2).
Target Name	APOH
Other Names	APOH;B2G1;B2GP1;BG
Accession No.	Uniprot:P02749GeneID:350
Uniprot	P02749
GeneID	350
SDS-PAGE MW	45kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Western blot analysis of extracts of 22Rv1 cells, using APOH antibody.

Background

Apolipoprotein H has been implicated in a variety of physiologic pathways including lipoprotein metabolism, coagulation, and the production of antiphospholipid autoantibodies. APOH may be a required cofactor for anionic phospholipid binding by the antiphospholipid autoantibodies found in sera of many patients with lupus and primary antiphospholipid syndrome, but it does not seem to be required for the reactivity of antiphospholipid autoantibodies associated with infections.

Note: This product is for in vitro research use only