

MAPK9 Antibody

Catalog No: #32258

Package Size: #32258-1 50ul #32258-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MAPK9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total MAPK9 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human JNK2 (NP_002743.3).
Target Name	MAPK9
Other Names	MAPK9;JNK-55;JNK2;JNK2A;JNK2ALPHA;JNK2B;JNK2BETA;PRKM9;SAPK;SAPK1a;p54a;p54aSAPK
Accession No.	Uniprot:P45984GeneID:5601
Uniprot	P45984
GeneID	5601
SDS-PAGE MW	54kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using JNK2 antibody.

Background

The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase targets specific transcription factors, and thus mediates immediate-early gene expression in response to various cell stimuli. It is most closely related to MAPK8, both of which are involved in UV radiation induced apoptosis, thought to be related to the cytochrome c-mediated cell death pathway. This gene and MAPK8 are also known as c-Jun N-terminal kinases. This kinase blocks the ubiquitination of tumor suppressor p53, and thus it increases the stability of p53 in nonstressed cells. Studies of this gene's mouse counterpart suggest a key role in T-cell differentiation. Several alternatively spliced transcript variants encoding distinct isoforms have been reported.

Note: This product is for in vitro research use only