

BrdU Mouse Monoclonal Antibody HRP Conjugated(A7)

Catalog No: #C08394H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	BrdU Mouse Monoclonal Antibody HRP Conjugated(A7)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	A7
Isotype	IgG
Purification	Purified by Protein G.
Applications	IHC-P IHC-F ICC
Species Reactivity	Rt
Crossing Reactivity	BrdU
Immunogen Description	BrdU conjugated to KLH
Conjugates	HRP
Target Name	BrdU
Other Names	Bromodeoxyuridine; Bromodeoxyuridine; 5-Bromo-2a2-deoxyuridine; 5-BrdU; proliferation Marker.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Background

Bromodeoxyuridine is a synthetic nucleoside that is an analogue of thymidine. BrdU is incorporated into DNA during S-phase of a cell cycle. It is commonly used in the detection of proliferating cells in living tissues

Note: This product is for in vitro research use only