

# CEA Mouse Monoclonal Antibody FITC Conjugated(B5 )

Catalog No: #C08406F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CEA Mouse Monoclonal Antibody FITC Conjugated(B5 )
Host Species	Mouse
Clonality	Monoclonal
Clone No.	B5
Isotype	IgG
Purification	Purified by Protein G.
Applications	IF
Species Reactivity	Hu
Immunogen Description	Human Carcinoembryonic Antigen
Conjugates	FITC
Target Name	CEA
Other Names	CEA; CD66e; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; Meconium antigen 1; CEACAM5
Accession No.	Swiss-Prot#P06731NCBI Gene ID1048
Uniprot	P06731
GeneID	1048;
Excitation Emission	494nm 518nm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins.

Note: This product is for in vitro research use only