

CEA Mouse Monoclonal Antibody HRP Conjugated(B5)

Catalog No: #C08406H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CEA Mouse Monoclonal Antibody HRP Conjugated(B5)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	B5
Isotype	IgG
Purification	Purified by Protein G.
Applications	IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	Human Carcinoembryonic Antigen
Conjugates	HRP
Target Name	CEA
Other Names	CEA; CD66e; Carcinoembryonic antigen-related cell adhesion molecule 5; Carcinoembryonic antigen; Meconium antigen 1; CEACAM5
Accession No.	Swiss-Prot#P06731NCBI Gene ID1048
Uniprot	P06731
GeneID	1048;
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Cell surface glycoprotein that plays a role in cell adhesion and in intracellular signaling. Receptor for E.coli Dr adhesins.

Note: This product is for in vitro research use only