

# RBC Mouse Monoclonal Antibody HRP Conjugated(G2H4 )

Catalog No: #C08414H

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RBC Mouse Monoclonal Antibody HRP Conjugated(G2H4 )
Host Species	Mouse
Clonality	Monoclonal
Clone No.	G2H4
Isotype	IgG
Purification	Purified by Protein G.
Applications	IHC-P IHC-F
Species Reactivity	Hu
Immunogen Description	Human washed pooled human Red Blood Cell
Conjugates	HRP
Target Name	RBC
Other Names	red blood cell; Erythrocytes; RBCs.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200

## Background

Red blood cells are the most common type of blood cell and are the vertebrate body's principal means of delivering oxygen from the lungs or gills to body tissues via the blood.

Note: This product is for in vitro research use only