

CSNK2A2 Antibody

Catalog No: #32343

Package Size: #32343-1 50ul #32343-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CSNK2A2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CSNK2A2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human CSNK2A2 (NP_001887.1).
Target Name	CSNK2A2
Other Names	CSNK2A2;CK2A2;CK2alpha;CSNK2A1
Accession No.	Uniprot:P19784GeneID:1459
Uniprot	P19784
GeneID	1459
SDS-PAGE MW	40kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200


Images




Western blot analysis of extracts of HeLa cells, using CSNK2A2 antibody.



Western blot analysis of extracts of various cell lines, using CSNK2A2 antibody.



Immunofluorescence analysis of PC-12 cells using CSNK2A2 Rabbit pAb.



Immunofluorescence analysis of NIH/3T3 cells using CSNK2A2 Rabbit pAb.

Note: This product is for in vitro research use only