

SFRS1 Antibody

Catalog No: #32366

Package Size: #32366-1 50ul #32366-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SFRS1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total SFRS1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human SRSF1/SF2/ASF (NP_008855.1).
Target Name	SFRS1
Other Names	SRSF1;ASF;SF2;SF2p33;SFRS1;SRp30a
Accession No.	Uniprot:Q07955GeneID:6426
Uniprot	Q07955
GeneID	6426
SDS-PAGE MW	21KDa/28KDa/33KDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SRSF1/SF2/ASF antibody.

Background

This gene encodes a member of the arginine/serine-rich splicing factor protein family. The encoded protein can either activate or repress splicing, depending on its phosphorylation state and its interaction partners. Multiple transcript variants have been found for this gene. There is a pseudogene of this gene on chromosome 13.

Note: This product is for in vitro research use only