

CXADR Antibody

Catalog No: #32451

Package Size: #32451-1 50ul #32451-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CXADR Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CXADR protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human CXADR.
Target Name	CXADR
Other Names	CAR; HCAR; CAR4/6;
Accession No.	Swiss-Prot:P78310NCBI Gene ID:1525
Uniprot	P78310
GeneID	1525;
SDS-PAGE MW	40KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CXADR at 1:1000 dilution.

Background

The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21.

Note: This product is for in vitro research use only