

FANCD2 Antibody

Catalog No: #32584

Package Size: #32584-1 50ul #32584-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FANCD2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Rat
Specificity	The antibody detects endogenous level of total FANCD2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human FANCD2 (NP_149075.2).
Target Name	FANCD2
Other Names	FA4;FAD;FACD;FAD2;FA-D2;FANCD;FANCD2
Accession No.	Uniprot:Q9BXW9GenelD:2177
Uniprot	Q9BXW9
GenelD	2177
SDS-PAGE MW	180kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IF □ 1:50 - 1:200

Images



Immunofluorescence analysis of C6 cells using FANCD2 Rabbit pAb.

