

## ITGB5 Antibody

Catalog No: #32674

Package Size: #32674-1 50ul #32674-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

|                       |  |
|-----------------------|--|
| Product Name          | ITGB5 Antibody   |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Purification          | Antibodies were purified by affinity purification using immunogen.   |
| Applications          | WB   |
| Species Reactivity    | Human,Mouse,Rat  |
| Specificity           | The antibody detects endogenous level of total ITGB5 protein.  |
| Immunogen Type        | Recombinant Protein  |
| Immunogen Description | Recombinant protein of human ITGB5.  |
| Target Name           | ITGB5  |
| Accession No.         | Swiss-Prot:P18084NCBI Gene ID:3693   |
| Uniprot               | P18084   |
| GeneID                | 3693;  |
| SDS-PAGE MW           | 88KD   |
| Concentration         | 1.0mg/ml   |
| Formulation           | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage               | Store at -20°C   |

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using ITGB5 at 1:500 dilution.

## Background

Integrin alpha-V/beta-5 is a receptor for fibronectin. It recognizes the sequence R-G-D in its ligand.

---

Note: This product is for in vitro research use only