

## SNCG Antibody

Catalog No: #32690

Package Size: #32690-1 50ul #32690-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	SNCG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total SNCG protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human gamma-Synuclein (NP_003078.2).
Target Name	SNCG
Other Names	SNCG;BCSG1;SR
Accession No.	Uniprot:O76070GeneID:6623
Uniprot	O76070
GeneID	6623
SDS-PAGE MW	13kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using gamma-Synuclein antibody.

## Background

---

This gene encodes a member of the synuclein family of proteins which are believed to be involved in the pathogenesis of neurodegenerative diseases. Mutations in this gene have also been associated with breast tumor development.

---

Note: This product is for in vitro research use only