

CFLAR Antibody

Catalog No: #32710

Package Size: #32710-1 50ul #32710-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CFLAR Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total CFLAR protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human cFLIP (NP_001120655.1).
Target Name	CFLAR
Other Names	CFLAR;CASH;CASP8AP1;CLARP;Casper;FLAME;FLAME-1;FLAME1;FLIP;I-FLICE;MRIT;c-FLIP;c-FLIPL;c-F LIPR;c-FLIPS
Accession No.	Uniprot:O15519GeneID:8837
Uniprot	O15519
GeneID	8837
SDS-PAGE MW	28kDa/55kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using cFLIP antibody.



Western blot analysis of extracts of Jurkat cells, using cFLIP antibody.

Background

The protein encoded by this gene is a regulator of apoptosis and is structurally similar to caspase-8. However, the encoded protein lacks caspase activity and appears to be itself cleaved into two peptides by caspase-8. Several transcript variants encoding different isoforms have been found for this gene, and partial evidence for several more variants exists.

Note: This product is for in vitro research use only