

## GNB2L1 Antibody

Catalog No: #32713

Package Size: #32713-1 50ul #32713-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	GNB2L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total GNB2L1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human RACK1 (NP_006089.1).
Target Name	GNB2L1
Other Names	RACK1;GNB2L1;Gnb2-rs1;H12.3;HLC-7;PIG21
Accession No.	Uniprot:P63244GeneID:10399
Uniprot	P63244
GeneID	10399
SDS-PAGE MW	35kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using RACK1 antibody.

---

Note: This product is for in vitro research use only