

BTC Antibody

Catalog No: #32736

Package Size: #32736-1 50ul #32736-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	BTC Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total BTC protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human BTC.
Target Name	BTC
Accession No.	Swiss-Prot:P35070NCBI Gene ID:685
Uniprot	P35070
GeneID	685;
SDS-PAGE MW	20KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using BTC at 1:1000 dilution.

Background

The protein encoded by this gene is a member of the EGF family of growth factors. It is synthesized primarily as a transmembrane precursor, which is

then processed to mature molecule by proteolytic events. This protein is a ligand for the EGF receptor.

Note: This product is for in vitro research use only