

## PTGES3 Antibody

Catalog No: #32773

Package Size: #32773-1 50ul #32773-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PTGES3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total PTGES3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human PTGES3.
Target Name	PTGES3
Other Names	P23; TEBP; cPGES;
Accession No.	Swiss-Prot:Q15185NCBI Gene ID:10728
Uniprot	Q15185
GeneID	10728;
SDS-PAGE MW	19KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using PTGES3 at 1:1000 dilution.

## Background

---

Molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes.

---

Note: This product is for in vitro research use only