

## PVRL2 Antibody

Catalog No: #32816

Package Size: #32816-1 50ul #32816-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PVRL2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total PVRL2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human Nectin 2/CD112 (NP_002847.1).
Target Name	PVRL2
Other Names	NECTIN2;CD112;HVEB;PRR2;PVRL2;PVRR2;nectin-2
Accession No.	Uniprot:Q92692GeneID:5819
Uniprot	Q92692
GeneID	5819
SDS-PAGE MW	72kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

## Images



Western blot analysis of extracts of various cell lines, using Nectin 2/CD112 antibody.

## Background

---

This gene encodes a single-pass type I membrane glycoprotein with two Ig-like C2-type domains and an Ig-like V-type domain. This protein is one of the plasma membrane components of adherens junctions. It also serves as an entry for certain mutant strains of herpes simplex virus and pseudorabies virus, and it is involved in cell to cell spreading of these viruses. Variations in this gene have been associated with differences in the severity of multiple sclerosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

---

Note: This product is for in vitro research use only