

## RANGAP1 Antibody

Catalog No: #32818

Package Size: #32818-1 50ul #32818-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	RANGAP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Specificity	The antibody detects endogenous level of total RANGAP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human RanGAP1 (NP_002874.1).
Target Name	RANGAP1
Other Names	RANGAP1;Fug1;RANGAP;SD
Accession No.	Uniprot:P46060GeneID:5905
Uniprot	P46060
GeneID	5905
SDS-PAGE MW	64kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

## Application Details

WB □ 1:500 - 1:2000

## Images



Western blot analysis of extracts of various cell lines, using RanGAP1 antibody.

## Background

---

This gene encodes a protein that associates with the nuclear pore complex and participates in the regulation of nuclear transport. The encoded protein interacts with Ras-related nuclear protein 1 (RAN) and regulates guanosine triphosphate (GTP)-binding and exchange. Alternative splicing results in multiple transcript variants.

---

Note: This product is for in vitro research use only