

PRKRA Antibody

Catalog No: #32843

Package Size: #32843-1 50ul #32843-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PRKRA Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse
Specificity	The antibody detects endogenous level of total PRKRA protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human PRKRA/PACT (NP_003681.1).
Target Name	PRKRA
Other Names	PRKRA;DYT16;HSD14;PACT;RAX
Accession No.	Uniprot:O75569GeneID:8575
Uniprot	O75569
GeneID	8575
SDS-PAGE MW	35kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200

Images



Western blot analysis of extracts of mouse liver, using PRKRA/PACT antibody.

Background

This gene encodes a protein kinase activated by double-stranded RNA which mediates the effects of interferon in response to viral infection. Mutations in this gene have been associated with dystonia. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only