

SNRPE Antibody

Catalog No: #32883

Package Size: #32883-1 50ul #32883-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SNRPE Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total SNRPE protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human SNRPE (NP_003085.1).
Target Name	SNRPE
Other Names	SNRPE;HYPT11;SME;Sm-E;snRNP-E
Accession No.	Uniprot:P62304GenelD:6635
Uniprot	P62304
GenelD	6635
SDS-PAGE MW	11kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded rat ovary using SNRPE Rabbit pAb.



Western blot analysis of extracts of various cell lines, using SNRPE antibody.



Immunohistochemistry of paraffin-embedded human thyroid cancer using SNRPE Rabbit pAb.

Background

The protein encoded by this gene is a core component of U small nuclear ribonucleoproteins, which are key components of the pre-mRNA processing spliceosome. The encoded protein plays a role in the 3' end processing of histone transcripts. This protein is one of the targets in the autoimmune disease systemic lupus erythematosus, and mutations in this gene have been associated with hypotrichosis. Several pseudogenes of this gene have been identified.

Note: This product is for in vitro research use only