

TLE1 Antibody

Catalog No: #32891

Package Size: #32891-1 50ul #32891-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TLE1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Specificity	The antibody detects endogenous level of total TLE1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant fusion protein of human TLE1 (NP_005068.2).
Target Name	TLE1
Other Names	TLE1;ESG;ESG1;GRG1
Accession No.	Uniprot:Q04724GeneID:7088
Uniprot	Q04724
GeneID	7088
SDS-PAGE MW	83kDa
Concentration	1.0mg/ml
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of THP-1 cells, using TLE1 antibody.

Background

Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. Enhances FOXG1/BF-1- and HES1-mediated transcriptional repression (By similarity. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG.

Note: This product is for in vitro research use only