

## GABARAP Antibody

Catalog No: #32924

Package Size: #32924-1 50ul #32924-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GABARAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Specificity	The antibody detects endogenous level of total GABARAP protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human GABARAP.
Target Name	GABARAP
Other Names	MM46; ATG8A; GABARAP-a;
Accession No.	Swiss-Prot:O95166NCBI Gene ID:11337
Uniprot	O95166
GeneID	11337;
SDS-PAGE MW	13KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


## Application Details

WB □ 1:500 - 1:2000 IHC □ 1:50 - 1:200 IF □ 1:50 - 1:200


## Images




Immunohistochemistry of paraffin-embedded rat kidney using GABARAP at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded mouse kidney using GABARAP at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat brain using GABARAP at dilution of 1:100 (40x lens).



Immunofluorescence analysis of MCF-7 cells using GABARAP . Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using GABARAP at 1:1000 dilution.

## Background

Gamma-aminobutyric acid A receptors.

Note: This product is for in vitro research use only