

ETV4 Antibody

Catalog No: #33049

Package Size: #33049-1 50ul #33049-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ETV4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were purified by affinity purification using immunogen.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total ETV4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant protein of human ETV4.
Target Name	ETV4
Other Names	E1AF; PEA3; E1A-F; PEAS3;
Accession No.	Swiss-Prot:P43268NCBI Gene ID:2118
Uniprot	P43268
GeneID	2118;
SDS-PAGE MW	53KD
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

Western blotting: 1:500 - 1:2000


Immunohistochemistry: 1:50 - 1:200

Immunofluorescence: 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using Pea3 / ETV4 antibody at 1:1000 dilution.



Immunofluorescence analysis of U2OS cells using Pea3 / ETV4 antibody. Blue: DAPI for nuclear staining.

Background

Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3'.

Note: This product is for in vitro research use only