

ABHD14B Antibody

Catalog No: #34369

Package Size: #34369-1 50ul #34369-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD14B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ABHD14B protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human ABHD14B.
Target Name	ABHD14B
Other Names	Abhydrolase domain-containing protein 14B; EC 3.-.-.-; CCG1-interacting factor B;
Accession No.	Swiss-Prot: Q96IU4NCBI Gene ID: 84836
Uniprot	Q96IU4
GeneID	84836;
SDS-PAGE MW	22kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue using ABHD14B antibody #34369.

Background

Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

Zilfou JT, et al. (2001) Mol Cell Biol 21(12): 3974-3985.

Boehm M, et al. (2002) EMBO J 21(13): 3390-3401.

Vadlamudi RK, et al. (2005) Mol Cell Biol 25(9): 3726-3736.

Note: This product is for in vitro research use only