

## BCLAF1 Antibody

Catalog No: #34481

Package Size: #34481-1 50ul #34481-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BCLAF1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total BCLAF1 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human BCLAF1.
Target Name	BCLAF1
Other Names	Bcl-2-associated transcription factor 1; BCLF1; Btf; KIAA0164;
Accession No.	Swiss-Prot: Q9NYF8NCBI Gene ID: 9774
Uniprot	Q9NYF8
GeneID	9774;
SDS-PAGE MW	106kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:500~1:3000

## Images



Western blot analysis of extracts from 293 cells, Jurkat cells and A549 cells, using BCLAF1 antibody #34481.

## Background

---

Death-promoting transcriptional repressor. May be involved in cyclin-D1/CCND1 mRNA stability through the SNARP complex which associates with both the 3'end of the CCND1 gene and its mRNA.

Kasof G.M., Mol. Cell. Biol. 19:4390-4404(1999).

Nagase T., DNA Res. 3:17-24(1996).

Mungall A.J., Nature 425:805-811(2003)

---

Note: This product is for in vitro research use only