

CBLN3 Antibody

Catalog No: #34572

Package Size: #34572-1 50ul #34572-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CBLN3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total CBLN3 protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human CBLN3.
Target Name	CBLN3
Other Names	Cerebellin-3;
Accession No.	Swiss-Prot: Q6UW01NCBI Gene ID: 643866
Uniprot	Q6UW01
GeneID	643866;
SDS-PAGE MW	22kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50~1:100

Images



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using CBLN3 antibody #34572.

Background

May be involved in synaptic functions in the CNS By similarity.

Clark H.F., Genome Res. 13:2265-2270(2003).

Zhang Z., Protein Sci. 13:2819-2824(2004).

Note: This product is for in vitro research use only