

SFRS2B Antibody

Catalog No: #35075



Package Size: #35075-1 50ul #35075-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SFRS2B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SFRS2B protein.
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from internal of human SFRS2B.
Target Name	SFRS2B
Other Names	SFR2B; SRp46 splicing factor;
Accession No.	Swiss-Prot: Q9BRL6NCBI Gene ID: 10929
Uniprot	Q9BRL6
GeneID	10929;
SDS-PAGE MW	32kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C


Application Details

WB 1:500-1:2000

IHC 1:100-1:300

IF 1:50-200

Images



Western blot analysis of extracts from HeLa cells, using SFRS2B antibody #35075.

Background

Involved in pre-mRNA alternative splicing.

Soret J., Mol. Cell. Biol. 18:4924-4934(1998).

Mural R.J., Submitted (JUL-2005) to the EMBL/GenBank/DDBJ databases.

Olsen J.V., Cell 127:635-648(2006).

Note: This product is for in vitro research use only