

## NEGR1 antibody

Catalog No: #22049

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	NEGR1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 82 and 341 of Human NEGR1
Target Name	NEGR1
Accession No.	Swiss-Prot:Q7Z3B1Gene ID:257194
Uniprot	Q7Z3B1
GeneID	257194;
Concentration	0.3mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 39kd

Western blotting: 1:500-1:3000

## Images



Sample(30 ug whole cell lysate)  
A: HeLa S3  
10% SDS PAGE  
Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only