

TAX1BP1 Antibody

Catalog No: #35972

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAX1BP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total TAX1BP1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human Tax1 (human T-cell leukemia virus type I) binding protein 1
Target Name	TAX1BP1
Other Names	T6BP; TXBP151; CALCOCO3
Accession No.	Swiss-Prot#: Q86VP1NCBI Gene ID: 8887Gene Accssion: BC050358
Uniprot	Q86VP1
GeneID	8887;
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:300

Images



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using #35972 at dilution 1/60.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #35972 at dilution 1/60.

Background

This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[]

Note: This product is for in vitro research use only