

SCGB1A1 Antibody

Catalog No: #35982

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SCGB1A1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SCGB1A1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human secretoglobin, family 1A, member 1 (uteroglobin)
Target Name	SCGB1A1
Other Names	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
Accession No.	Swiss-Prot#: P11684NCBI Gene ID: 7356Gene Accssion: BC004481
Uniprot	P11684
GeneID	7356;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #35982 at dilution 1/60.

Background

This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a

susceptibility to asthma.

Note: This product is for in vitro research use only