

## ADCK4 Antibody

Catalog No: #36049

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ADCK4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ADCK4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human aarF domain containing kinase 4
Target Name	ADCK4
Other Names	FLJ12229;aarF domain containing kinase 4
Accession No.	Swiss-Prot#: Q96D53NCBI Gene ID: 79934Gene Accssion: BC013114
Uniprot	Q96D53
GeneID	79934;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:25-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36049 at dilution 1/30.

## Background

This gene encodes a protein with two copies of a domain found in protein kinases. The encoded protein has a complete protein kinase catalytic domain, and a truncated domain that contains only the active and binding sites of the protein kinase domain, however, it is not known whether the

protein has any kinase activity. Multiple transcript variants encoding different isoforms have been found for this gene.?

---

Note: This product is for in vitro research use only