

AGPAT4 Antibody

Catalog No: #36063

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AGPAT4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AGPAT4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human 1-acylglycerol-3-phosphate O-acyltransferase 4
Target Name	AGPAT4
Other Names	1-AGPAT4; dJ473J16.2; LPAAT-delta; RP3-473J16.2
Accession No.	Swiss-Prot#: Q9NRZ5NCBI Gene ID: 56895Gene Accssion: BC020209
Uniprot	Q9NRZ5
GeneID	56895;
Concentration	1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #36063 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36063 at dilution 1/25.

Background

This gene encodes a member of the 1-acylglycerol-3-phosphate O-acyltransferase family. This integral membrane protein converts lysophosphatidic acid to phosphatidic acid, the second step in de novo phospholipid biosynthesis.

Note: This product is for in vitro research use only