

## AGTRAP Antibody

Catalog No: #36065

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AGTRAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total AGTRAP protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human angiotensin II receptor-associated protein
Target Name	AGTRAP
Other Names	ATRAP
Accession No.	Swiss-Prot#: Q6RW13NCBI Gene ID: 57085Gene Accsion: BC017328
Uniprot	Q6RW13
GeneID	57085;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36065 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #36065 at dilution 1/40.

## Background

This gene encodes a transmembrane protein localized to the plasma membrane and perinuclear vesicular structures. The gene product interacts with the angiotensin II type I receptor and negatively regulates angiotensin II signaling. Alternative splicing of this gene generates multiple transcript variants encoding different isoforms.

Note: This product is for in vitro research use only