

ALG9 Antibody

Catalog No: #36081

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ALG9 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ALG9 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human ALG9, alpha-1,2-mannosyltransferase
Target Name	ALG9
Other Names	CDG1L; DIBD1; LOH11CR1J
Accession No.	Swiss-Prot#: Q9H6U8NCBI Gene ID: 79796Gene Accssion: BC009255
Uniprot	Q9H6U8
GeneID	79796;
Concentration	2.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #36081 at dilution 1/60.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36081 at dilution 1/60.

Background

This gene encodes an alpha-1,2-mannosyltransferase enzyme that functions in lipid-linked oligosaccharide assembly. Mutations in this gene result in congenital disorder of glycosylation type II. Multiple transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only