

APOBEC4 Antibody

Catalog No: #36120

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	APOBEC4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total APOBEC4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 4 (putative)
Target Name	APOBEC4
Other Names	C1orf169
Accession No.	Swiss-Prot#: Q8WW27NCBI Gene ID: 403314Gene Accssion: BC021711
Uniprot	Q8WW27
GeneID	403314;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #36120 at dilution 1/20.

Background

This gene encodes a member of the AID/APOBEC family of polynucleotide (deoxy)cytidine deaminases, which convert cytidine to uridine. Other AID/APOBEC family members are involved in mRNA editing, somatic hypermutation and recombination of immunoglobulin genes, and innate

immunity to retroviral infection.?

Note: This product is for in vitro research use only