

ARL11 Antibody

Catalog No: #36261

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ARL11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ARL11 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	ARL11
Other Names	ARLTS1
Accession No.	Swiss-Prot#: Q969Q4NCBI Gene ID: 115761Gene Accssion: BC013150
Uniprot	Q969Q4
GeneID	115761;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36261 at dilution 1/30.

Background

This gene encodes a tumor suppressor related to the ADP-ribosylation factor (ARF) family of proteins. The encoded protein may play a role in apoptosis in a caspase-dependent manner. Polymorphisms in this gene have been associated with some familial cancers.

Note: This product is for in vitro research use only