

## ASB3 Antibody

Catalog No: #36264

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	ASB3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ASB3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human ankyrin repeat and SOCS box containing 3
Target Name	ASB3
Other Names	ASB-3
Accession No.	Swiss-Prot#: Q9Y575NCBI Gene ID: 51130Gene Accssion: BC006488/Q9Y575
Uniprot	Q9Y575
GeneID	100302652;51130;
SDS-PAGE MW	58kd
Concentration	0.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

## Images



Gel: 8%SDS-PAGE  
 Lysate: 40ug A549 cell  
 Primary antibody: 1/200 dilution  
 Secondary antibody dilution: 1/8000  
 Exposure time: 10 seconds

## Background

The protein encoded by this gene is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat

sequence and SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. Alternatively spliced transcript variants have been described for this gene.?

---

Note: This product is for in vitro research use only