

PGCP antibody

Catalog No: #22111

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	PGCP antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 19 and 266 of Human PGCP
Target Name	PGCP
Accession No.	Swiss-Prot:Q9Y646Gene ID:10404
Uniprot	Q9Y646
GeneID	10404;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 52kd

Western blotting: 1:500-1:3000

Images



Sample(30 ug of whole cell lysate)

A: Raji

7.5% SDS PAGE

Primary antibody diluted at 1: 500

Background

Carboxypeptidase that may play an important role in the hydrolysis of circulating peptides.

Note: This product is for in vitro research use only