

CHMP4B Antibody

Catalog No: #36350

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHMP4B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total CHMP4B protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human charged multivesicular body protein 4B
Target Name	CHMP4B
Other Names	SNF7; CTPP3; Shax1; CHMP4A; SNF7-2; VPS32B; CTRCT31; Vps32-2; C20orf178; dJ553F4.4
Accession No.	Swiss-Prot#: Q9H444NCBI Gene ID: 128866Gene Accssion: BC033859
Uniprot	Q9H444
GeneID	128866;
Concentration	1.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36350 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #36350 at dilution 1/40.

Background

This gene encodes a member of the chromatin-modifying protein/charged multivesicular body protein (CHMP) protein family. The protein is part of the endosomal sorting complex required for transport (ESCRT) complex III (ESCRT-III), which functions in the sorting of endocytosed cell-surface receptors into multivesicular endosomes. The ESCRT machinery also functions in the final abscission stage of cytokinesis and in the budding of enveloped viruses such as HIV-1.

Note: This product is for in vitro research use only