

DLL4 Antibody

Catalog No: #36418

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DLL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DLL4 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human delta-like 4 (Drosophila)
Target Name	DLL4
Other Names	hdelta2
Accession No.	Swiss-Prot#: Q9NR61NCBI Gene ID: 54567Gene Accssion: BC106950
Uniprot	Q9NR61
GeneID	54567;
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #36418 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36418 at dilution 1/40.

Background

This gene is a homolog of the *Drosophila* delta gene. The delta gene family encodes Notch ligands that are characterized by a DSL domain, EGF repeats, and a transmembrane domain.

Note: This product is for in vitro research use only