

## FBP2 Antibody

Catalog No: #36472

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	FBP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FBP2 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	FBP2
Other Names	F16P2;FBPase 2;D-fructose-1
Accession No.	Swiss-Prot#: O00757NCBI Gene ID: 8789Gene Accssion: BC113632
Uniprot	O00757
GeneID	8789;
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #36472 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #36472 at dilution 1/40.

## Background

This gene encodes a gluconeogenesis regulatory enzyme which catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate.

Note: This product is for in vitro research use only