

GTPBP10 Antibody

Catalog No: #36516

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GTPBP10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GTPBP10 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human GTP-binding protein 10
Target Name	GTPBP10
Other Names	ObgH2; UG0751c10
Accession No.	Swiss-Prot#: A4D1E9NCBI Gene ID: 85865Gene Accssion: BC021573
Uniprot	A4D1E9
GeneID	85865;
SDS-PAGE MW	43kd
Concentration	1.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 8%SDS-PAGE
 Lysate: 40ug A549 cell
 Primary antibody: 1/200 dilution
 Secondary antibody dilution: 1/8000
 Exposure time: 10 minutes

Background

Small G proteins, such as GTPBP10, act as molecular switches that play crucial roles in the regulation of fundamental cellular processes such as protein synthesis, nuclear transport, membrane trafficking, and signal transduction

Note: This product is for in vitro research use only