

LCMT1 Antibody

Catalog No: #36582

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LCMT1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total LCMT1 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to residues near the C terminal of human leucine carboxyl methyltransferase 1
Target Name	LCMT1
Other Names	LCMT; PPMT1; CGI-68
Accession No.	Swiss-Prot#: Q9UIC8NCBI Gene ID: 51451Gene Accssion: BC014217
Uniprot	Q9UIC8
GeneID	51451;
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using #36582 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #36582 at dilution 1/20.

Background

LCMT1 catalyzes the methylation of the carboxyl group of the C-terminal leucine residue (leu309) of the catalytic subunit of protein phosphatase-2A.

Note: This product is for in vitro research use only