

## MAGEB10 Antibody

Catalog No: #36596

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	MAGEB10 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MAGEB10 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human melanoma antigen family B, 10
Target Name	MAGEB10
Other Names	FLJ32965;MAGBA;melanoma antigen family B; 10
Accession No.	Swiss-Prot#: Q96LZ2 NCBI Gene ID: 139422Gene Accssion: BC101308
Uniprot	Q96LZ2
GeneID	139422;
Concentration	2.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN3, 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:300

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36596 at dilution 1/60.



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #36596 at dilution 1/60.

## Background

This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells.

Note: This product is for in vitro research use only