

MBNL3 Antibody

Catalog No: #36603

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MBNL3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MBNL3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from internal residues of human muscleblind-like splicing regulator 3
Target Name	MBNL3
Other Names	CHCR; MBLX; MBXL; MBLX39
Accession No.	Swiss-Prot#: Q9NUK0NCBI Gene ID: 55796Gene Accssion: BC042090/Q9NUK0
Uniprot	Q9NUK0
GeneID	55796;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #36603 at dilution 1/30.



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #36603 at dilution 1/30.

Background

This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described.

Note: This product is for in vitro research use only