

cofilin(Phospho-Ser3) Antibody

Catalog No: #11139

Package Size: #11139-1 50ul #11139-2 100ul

Description

Product Name	cofilin(Phospho-Ser3) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB IHC IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of cofilin only when phosphorylated at serine 3.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of serine 3 (M-A-S(p)-G-V) derived from Human cofilin.
Conjugates	Unconjugated
Target Name	cofilin
Modification	Phospho
Other Names	CFL; CFL1;
Accession No.	Swiss-Prot: P23528NCBI Protein: NP_005498.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details


Predicted MW: 19kd

Western blotting: 1:500~1:1000


Immunohistochemistry: 1:50~1:100

Immunofluorescence: 1:100~1:200


Images




Western blot analysis of extracts from COLO205 cells using cofilin(Phospho-Ser3) Antibody #11139(Lane 2) and the same antibody preincubated with blocking peptide(Lane1).



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using cofilin(Phospho-Ser3) Antibody #11139(left) or the same antibody preincubated with blocking peptide(right).



Immunofluorescence staining of methanol-fixed HeLa cells using cofilin(Phospho-Ser3) Antibody #11139.



Lane 1 : HUVEC; Lane 2 : HUVEC treated with EGF; Lane 3 : MCF; Lane 4 : MCF treated with EGF; Lysates/proteins at 40 µg per lane. Predicted band size : 19 kDa Observed band size : 19 kDa

Background

Cofilin reversibly actin polymerization and depolymerization in a pH-sensitive manner. It has the ability to bind G- and F-actin in a 1:1 ratio of cofilin to actin. It is the major component of intranuclear and cytoplasmic actin rods.

Kobayashi M, et al. (2006) EMBO J 25(4): 713-26.

Wang Y, et al. (2005) Biol Chem 280(13): 12683-9.

Smith-Beckerman DM, et al. (2005) Mol Cell Proteomics: 156-68.

Guillaume Huet¹, Eeva Kaisa Rajakyla, Tiina Viita¹ et al., Actin-regulated feedback loop based on Phactr4, PP1 and cofilin maintains the actin monomer pool., *Journal of Cell Science*, 126:497η— C507(2012)

[PMID:23203801](#)

Hui Guo, Yi Lv, Tao Tian et al., Downregulation of p57 accelerates the growth and invasion of hepatocellular carcinoma., *Carcinogenesis*, vol.32 no.12 pp.1897η— C1904(2011)

[PMID:22002319](#)

Lijun Zhang, Jun Luo, Ping Wan et al., Regulation of cofilin phosphorylation and asymmetry in collective cell migration during morphogenesis., *Development.*, 138:455-464(2011)

[PMID:21205790](#)

et al., Downregulation of p57 Accelerates the Growth and Invasion of Hepatocellular Carcinoma. In *Carcinogenesis* on 2011 Dec by Hui Guo , Yi Lv, et al..PMID: 22002319, , (2011)

[PMID:22002319](#)

et al., Regulation of cofilin phosphorylation and asymmetry in collective cell migration during morphogenesis. In *Development* on 2011 Feb by Lijun Zhang, Jun Luo, et al..PMID:

21205790, , (2011)

[PMID:21205790](#)

et al., Actin-regulated feedback loop based on Phactr4, PP1 and cofilin maintains the actin monomer pool. In *J Cell Sci* on 2013 Jan 15 by Tiina Viita, Kari-Pekka Skarp, et al..PMID: 23203801, , (2013)

[PMID:23203801](#)

et al., The role of cytoplasmic p57 in invasion of hepatocellular carcinoma. In *BMC Gastroenterol.* On 2015 Aug 15 by Hui Guo, Yi Li et al..PMID:26271467, , (2015)

[PMID:
26271467](#)

et al., Genome-wide RNAi Screen for Nuclear Actin Reveals a Network of Cofilin Regulators .In *J Cell Sci* on 2015 Jul 1 by Joseph Dopie, Eeva K Rajakylⁿ ² et al..PMID: 26021350, , (2015)

[PMID:26021350](#)

et al., Suppression of p21Rac Signaling and Increased Innate Immunity Mediate Remission in Crohni ζ Disease. In *sci Transl Med* on 2014 Apr 23 by Kaushal Parikh , Lu Zhou et al..PMID:24760188, , (2014)

[PMID:24760188](#)

et al., The NAV2 homolog Sickie regulates F-actin-mediated axonal growth inθ• ½Drosophilaθ• ½mushroom body neurons via the non-canonical Rac-Cofilin pathway. In *Development* on 2014 Dec by Takashi Abe, Daisuke Yamazaki et al..PMID:25411210, , (2014)

[PMID:25411210](#)

et al., X chromosome-linked intellectual disability protein PQBP1 associates with and regulates the translation of specific mRNAs. In *Hum Mol Genet* on 2015 Aug 15 by Didi Wan, Zi Chao Zhang, et al..PMID:26002102, , (2015)

[PMID:26002102](#)

et al., Scribble scaffolds a signalosome for active forgetting. In *Neuron* on 2016 Jun 15 by Isaac Cervantes-Sandoval, Molee Chakraborty, et al..PMID:27263975

, , (2016)

[PMID:27263975](#)

et al., LIM Kinase Regulation of Cytoskeletal Dynamics Is Required for Salivary Gland Branching Morphogenesis. In *Mol Biol Cell* on 2014 Aug 15 by Shayoni Ray, Joseph A Fantii ε t al..PMID:24966172, , (2014)

[PMID:24966172](#)

el at., Regulation of nuclear actin levels and MRTF/SRF target gene expression during PC6.3 cell differentiation. In Exp Cell Res on 2022 Nov 15 by Salla Kyher?inen, Alise Hyrskyluoto, et al..PMID: 36122768, , (2022)

[PMID:36122768](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.