

GAP43(Phospho-Ser41) Antibody

Catalog No: #11281

Package Size: #11281-1 50ul #11281-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	GAP43(Phospho-Ser41) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	IF
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of GAP43 only when phosphorylated at Ser41.
Immunogen Type	Peptide-KLH
Immunogen Description	Peptide sequence around phosphorylation site of Ser41 (Q-A-S(p)-F-R) derived from Human GAP43.
Conjugates	Unconjugated
Target Name	GAP43
Modification	Phospho
Other Names	B-50; BASP2; NEUM; PP46; axonal membrane protein GAP-43
Accession No.	Swiss-Prot: P17677NCBI Protein: NP_001123536.1
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 43kd

Immunofluorescence: 1:100~1:200

Images



Immunofluorescence staining of methanol-fixed HeLa cells using GAP43(Phospho-Ser41) Antibody #11281.

Background

GAP43 encoded by this gene has been termed a 'growth' or 'plasticity' protein because it is expressed at high levels in neuronal growth cones during development and axonal regeneration. This protein is considered a crucial component of an effective regenerative response in the nervous system. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Rachael L. Neve, et.al. (1998) J. Neurosci; 18: 7757.

Yiping Shen, et.al. (2002) J. Neurosci; 22: 239.

Chantal Gamby, et.al. (1996) J. Biol. Chem; 271: 26698.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.